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Urgent Security Information

2018

Possible risks related to electronic treatment tables.

Important information for therapy and treatment table users and operators.

Identification of affected products

All treatment tables of the following series with a MAX60Motor:

Mahe E, MAUI E, SAMOA E, SAMOA EE, BALI E, GOZO E1 – E3, Vojta Bobath E

Dear customer,
we would like to inform you about the possible risk of your electronic treatment table.

Description of the problem and determined cause

Electrically adjustable treatment tables which are older than 10 years may drop down uncontrollably.

In case a patient is lying on it at that moment, he or she could be startled or at worst, fall down and be injured.

The problem occurs due to the SKF electric MAX60 motors lifespan of 10 years or 10.000 double strokes.

It is possible that the lifting motors are impaired by the ageing of the mechanics, abrasion or the motors mounting. This may lead to a sudden drop of the framework.

This concerns all treatment tables with a MAX60 lifting motor manufactured by SKF (formerly Magnetic).

What measures are to be applied and how can you avoid possible risks?

As an inspection is not provided by the manufacturer and spare parts are not available, all lifting motors older than 10 years should be replaced. Even if there is no visible abrasion or damage.

Please contact us for replacement of the lifting motor.
We will gladly create an individual offer for you.

Until the replacement we recommend not to use the lifting motor of respective treatment tables and to leave the table in a low position.

The ultimate decision how to handle the risk of the drive is up to you.

Note

Please note that as an operator, you are responsible for the security of your medical devices according to § 4 Medical Device Law (MPG) and § 2 Medical Devices Operator Regulation (MPBetreibV).

[Transfer of the information described](#)

Please ensure within your organization that all users receive this urgent security information. If you have passed on products to a third party, please forward this information or notify us by sending an e-mail.

The Federal Institute of Drugs and Medical Devices has received a copy of this 'Urgent Security Information'.

[Contact](#)

For questions regarding this security information please contact:

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Physioline GmbH Med. Praxiseinrichtungen